

Dissolution Tester DP-C8

Product Description



Features:

1. Auto-diagnose, auto-alarm.
2. Stainless steel basket and paddle.
3. Micro processor and LCD Display.
4. The head part of the instrument can be turned over smoothly and flexibly (Rise and fall for DP-C8).
5. Water circulating in the bath to keep the temperature conformity with magnetic pump.

Model	DP-C8
Basket	8pcs
Speed Range	20~200rpm
Speed Precision	±2rpm
Temp. Range	RT +45.0°C
Temp. Precision	±0.3°C
Timing Range	1~999min
Time Precision	±0.5min
Working Time	>24h
Paddle Swing Amplitude	±0.5mm
Basket Swing Amplitude	±1.0mm
Power Consumption	1200W
Power Supply	AC110/220V±10%, 50/60Hz
External Size (W*D*H) mm	680*540*480
Package Size (W*D*H) mm	1010*680*580